

**Contact:**

Lori Melton, Director of Marketing and Communications

[lmelton@rhythmink.com](mailto:lmelton@rhythmink.com)

866.633.3754 (toll-free)

803.252.1222

803.252.1111 (fax)

**RHYTHMLINK INTERNATIONAL  
RELEASES SPIDER™ EEG ELECTRODE**

COLUMBIA, S.C. - June, 2008 - Rhythmink International, a company that designs, manufactures and distributes neurodiagnostic accessories and provides custom OEM and manufacturing solutions, recently released its innovative Spider EEG Electrode.

The Spider EEG Electrode provides Rhythmink's customers with an alternative to traditional cup electrodes and features a flat, web-like design that favors patient comfort and offers superior adhesion capabilities. Currently, Rhythmink offers its Spider Electrode in two sizes, 6 and 10 mm, and types, Gold and Silver/Silver-Chloride. All varieties are available in 1.0, 1.5, 2.0, and 2.5 meter leadwire lengths.

For more information about Rhythmink, its Spider EEG Electrode, other products or services, visit [Rhythmink.com](http://Rhythmink.com).

###

**About Rhythmink International, LLC**

Rhythmink International, LLC designs, manufactures and distributes neurodiagnostic accessories and provides custom OEM and manufacturing solutions to its customers. Based in Columbia, S.C., Rhythmink was established in 2002 and is owned and operated by CEO Shawn Regan and COO Michael O'Leary. Rhythmink's advancements and improvements in the areas of technology, business development and corporate branding have been recognized by entities such as the Columbia Information Technology Council, *Inc. Magazine* and the International Academy of the Visual Arts. For more information about Rhythmink International, its products or services, visit [Rhythmink.com](http://Rhythmink.com).

**Photographs**

Silver/Silver-Chloride Reusable Spider EEG Electrode in 10 millimeter.