

Contact:

Lori Melton, Director of Marketing and Communications
lmelton@rhythmink.com
866.633.3754 (toll-free)
803.252.1222
803.252.1111 (fax)



**TWO ARE BETTER THAN ONE:
RHYTHMLINK INTERNATIONAL INTRODUCES PAIRED SURFACE ELECTRODES**

COLUMBIA, S.C. - Rhythmink International, a company that designs, manufactures and distributes neurodiagnostic accessories and provides custom packaging, private labeling, custom products and contract manufacturing solutions, has recently added Paired Recording and Stimulating Sticky Pad™ Surface Electrodes to its commercial product offerings.

Equipped with two different specially-formulated hydrogels which minimize stimulation risks (Stimulating) and optimize signal recording (Recording), Rhythmink's Paired Stimulating and Recording Sticky Pad Electrodes are available in 15x20 and 20x25 millimeter pad sizes and 1.0, 1.5 and 2.5 meter leadwire lengths.

For more information about Rhythmink's Paired Sticky Pad Surface Electrodes, other products or services, visit www.Rhythmink.com.

###

About Rhythmink International, LLC

Rhythmink International, LLC designs, manufactures and distributes neurodiagnostic accessories and provides custom packaging, private labeling, custom products and contract manufacturing to its customers. Based in Columbia, S.C., Rhythmink was established in 2002 and is owned and operated by CEO Shawn Regan and COO Michael O'Leary. Rhythmink's advancements and improvements in the areas of technology, business development and corporate branding have been recognized by entities such as the Columbia Information Technology Council, *Inc. Magazine* and the International Academy of the Visual Arts. For more information about Rhythmink International, its products or services, visit Rhythmink.com.




Photo 1: Paired Stimulating (left) and Paired Recording (right) Sticky Pad Surface Electrodes provide exclusive features designed to product safe results during patient procedures.